



NeuroSpeed® PTA Balloon Catheter



Balloon marker
(Additional marker in the catheter tip)

REACH.
PASS.
RESTORE.

xcandis®

ENGINEERING STROKE SOLUTIONS


NeuroSpeed® PTA Balloon Catheter

COMPLIANCE CHART

	Inflation Pressure	Balloon Diameter (mm)					
	bar	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
	2.0	1.21	1.72	2.09	2.42	3.06	3.26
	4.0	1.37	1.84	2.33	2.78	3.25	3.72
Nominal Pressure	6.0	1.50	2.00	2.50	3.00	3.50	4.00
	8.0	1.67	2.16	2.65	3.22	3.69	4.23
	10.0	1.85	2.27	2.75	3.38	3.83	4.37
	12.0	2.02	2.39	2.87	3.54	3.97	4.53
Rated Burst Pressure*	14.0	2.20	2.52	2.98	3.73	4.11	4.70

* Do not exceed!

SPECIFICATIONS



ORDERING INFORMATION

Labelled NeuroSpeed® Dimensions (mm)	Reference Number	Balloon Diameter (mm)	Balloon Working Length (mm)	ID (Inch)	OD dist. / prox. (French)	Usable Length (cm)
1.5 × 8.0	01-000605	1.5	8.0	0.0165	2.7 / 3.7	150
2.0 × 8.0	01-000600	2.0	8.0	0.0165	2.7 / 3.7	150
2.5 × 8.0	01-000601	2.5	8.0	0.0165	2.7 / 3.7	150
3.0 × 8.0	01-000602	3.0	8.0	0.0165	2.7 / 3.7	150
3.5 × 8.0	01-000603	3.5	8.0	0.0165	2.7 / 3.7	150
4.0 × 8.0	01-000604	4.0	8.0	0.0165	2.7 / 3.7	150

All changes or modifications, may they be technical or other, or changes in the availability of products are expressively reserved.



CE 0297

ACANDIS GmbH & Co. KG
Theodor-Fahrner-Str. 6
75177 Pforzheim
Germany

Tel: +49 7231 155 00 0
Fax: +49 7231 155 00 129
E-Mail: info@acandis.com
www.acandis.com

Distributed by: